

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Mouse anti Human CD8 FITC - CD45RA PE, Clone: [17D8/BHKG-2], FITC/RPE, Monoclonal NMB-845RA

Artikelname	Mouse anti Human CD8 FITC - CD45RA PE, Clone: [17D8/BHKG-2], FITC/RPE, Monoclonal
Artikelnummer	NMB-845RA
Hersteller Artikelnummer	845RA
Alternativnummer	NMB-845RA
Hersteller	NordicMubio
Wirt	Mouse
Kategorie	Antikörper
Applikation	FC
Spezies Reaktivität	Human
Immunogen	CD45=Derived from the hybridization of mouse NS-1 myeloma cells with spleen cells of BALB/c mice immunized with peripheral blood monocytes from rheumatoid arthritis patient CD8=Derived from the hybridization of mouse NS-1 myeloma cells with spleen cells from BALB/c mice immunized with human perherial blood T lymphocytes.
Konjugation	FITC/RPE
Produktbeschreibung	CD8 is useful for the Identification of human T cells suppressor/cytotoxic expressing the 32 kDa M.W. surface antigen. While CD45RA recognizes MW. 220-kdalton of the leucocyte antigen (LCA). CD45RA antigen is a member of the CD45 antigen family that ...

Klonalität	Monoclonal
Konzentration	Titered for flow cytometry
Klon-Bezeichnung	[17D8/BHKG-2]
Isotyp	IgG1, IgG1
UniProt	P01732
Puffer	Provided as solution in phosphate buffered saline with 0.08% sodium azide and 0.2% carrier protein
Reinheit	Protein A/G Chromatography
Formulierung	Bi-Test (FITC/RPE) Reagent
Formel	Provided as solution in phosphate buffered saline with 0.08% sodium azide and 0.2% carrier protein
Anwendungsbeschreibung	Add 10 µl of MAB / 100 µl of 106 PBMC in PBS. Mix Wash twice with PBS and analyze or fix with 0.5% v/v of paraformaldehyde in PBS and analyze. Add 10 µl of MAB / 100 µl of Whole Blood. Mix gently and incubate for 15 minutes at room temperature 20oC. Lyse Whole Blood and analyze or fix with 0.5% v/v paraformaldehyde in PBS and analyze. See instrument manufacturers instructions for Lysed Whole Blood and Immunofluorescence analysis with a flow cytometer or microscope.