

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### **Melanoma Marker(HMB45 + M2-7C10 + M2-9E3 + T311), CF405S conjugate, 0.1mg/mL, Clone: [HMB45 M2-7C10 M2-9E3 T311], Mouse, Monoclonal BOT-BNC040701-100**

Artikelname	Melanoma Marker(HMB45 + M2-7C10 + M2-9E3 + T311), CF405S conjugate, 0.1mg/mL, Clone: [HMB45 M2-7C10 M2-9E3 T311], Mouse, Monoclonal
Artikelnummer	BOT-BNC040701-100
Hersteller Artikelnummer	BNC040701-100
Alternativnummer	BOT-BNC040701-100-100UL
Hersteller	Biotium
Wirt	Mouse
Kategorie	Antikörper
Applikation	IHC, WB
Spezies Reaktivität	Human
Immunogen	Recombinant hMART-1 protein (M2-7C10, M2-9E3), Recombinant tyrosinase protein (T311), Extract of pigmented melanoma metastases from lymph nodes (HMB45)
Konjugation	CF405S
Produktbeschreibung	This antibody cocktail recognizes three melanoma-specific proteins, which include MART-1, Tyrosinase and gp100. MART-1 is a newly identified melanocyte differentiation antigen recognized by autologous cytotoxic T lymphocytes. Tyrosinase is one of the...
Klonalität	Monoclonal

Konzentration	0.1 mg/mL
Klon-Bezeichnung	[HMB45 M2-7C10 M2-9E3 T311]
Molekulargewicht	20-22 kDa (doublet) (MART), 70-80 kDa (Tyrosinase), 90-100 kDa (gp100)
UniProt	<a href="#">Q16655</a>
Puffer	PBS, 0.1% BSA, 0.05% azide
Quelle	Animal
Anwendungsbeschreibung	Higher concentration may be required for direct detection using primary antibody conjugates than for indirect detection with secondary antibody Immunohistology (Formalin-fixed) 0.5-1.0 ug/mL Staining of formalin-fixed tissues is enhanced by boiling tissue sections in 10 mM citrate buffer, pH 6.0, for 10-20 min followed by cooling at RT for 20 minutes Optimal dilution for a specific application should be determined by user