

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### **SARS-CoV-2 (COVID-19) 3CL-Mpro Protein, Tag-free, E. coli TRZ-P2020-019\_1000**

|                          |                                                                                                                                                                                                                                                               |
|--------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Artikelname              | SARS-CoV-2 (COVID-19) 3CL-Mpro Protein, Tag-free, E. coli                                                                                                                                                                                                     |
| Artikelnummer            | TRZ-P2020-019_1000                                                                                                                                                                                                                                            |
| Hersteller Artikelnummer | P2020-019_1000                                                                                                                                                                                                                                                |
| Alternativnummer         | TRZ-P2020-019_1000                                                                                                                                                                                                                                            |
| Hersteller               | trenzyme                                                                                                                                                                                                                                                      |
| Wirt                     | E. coli                                                                                                                                                                                                                                                       |
| Kategorie                | Biochemikalien                                                                                                                                                                                                                                                |
| Applikation              | ELISA, WB                                                                                                                                                                                                                                                     |
| Spezies Reaktivität      | Virus                                                                                                                                                                                                                                                         |
| Produktbeschreibung      | The new coronavirus SARS-CoV-2 expresses two proteases, the papain-like protease (PLpro) and 3C-like protease (3CLpro). Both belong to the group of cysteine proteases, as they have a cysteine residue at their catalytic site. Their main function is th... |
| Molekulargewicht         | 36,0 kDa                                                                                                                                                                                                                                                      |
| UniProt                  | <a href="#">P0DTD2</a>                                                                                                                                                                                                                                        |
| Puffer                   | PBS, contains Glycerol as protectant                                                                                                                                                                                                                          |
| Reinheit                 | > 90% as determined by SDS-PAGE                                                                                                                                                                                                                               |
| Formulierung             | liquid                                                                                                                                                                                                                                                        |

|         |                                                                                                                                                                                                                                                                                          |
|---------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Sequenz | SGFRKMAFPSGKVEGCMVQVTCGTTTLNGLWLDDVVYCPRHVICTSEDMLN<br>PNYEDLLIR<br>KSNHNFLVQAGNVQLRVIGHSMQNCVLKLKVDATANPKTPKYKVVRIQPGQTF<br>SVLACYNG<br>SPSGVYQCAMRPNFTIKGSFLNGSCGSVGFNIDYDCVSFCYMHHMELPTGVH<br>AGTDLEGN<br>FYGPFVDRQTAQAAGTDTTITVNVLAWLYAAVINGDRWFLNRFTTTLNDFNL<br>VAMKYNYE PLTQDHVDIL |
| Formel  | pH 7,4                                                                                                                                                                                                                                                                                   |