

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### **GENLISA Human ATP-Sensitive Inward Rectifier Potassium Channel 14 (KCNJ14 / KIR2.4) ELISA KBS-KBH1810**

|                          |                                                                                                                                                                                                                                                               |
|--------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Artikelname              | GENLISA Human ATP-Sensitive Inward Rectifier Potassium Channel 14 (KCNJ14 / KIR2.4) ELISA                                                                                                                                                                     |
| Artikelnummer            | KBS-KBH1810                                                                                                                                                                                                                                                   |
| Hersteller Artikelnummer | KBH1810                                                                                                                                                                                                                                                       |
| Alternativnummer         | KBS-KBH1810-1                                                                                                                                                                                                                                                 |
| Hersteller               | Krishgen Biosystems                                                                                                                                                                                                                                           |
| Kategorie                | Kits/Assays                                                                                                                                                                                                                                                   |
| Spezies Reaktivität      | Human                                                                                                                                                                                                                                                         |
| Produktbeschreibung      | The GENLISA Human Atp-Sensitive Inward Rectifier Potassium Channel 14 (KCNJ14 / KIR2.4) ELISA includes features like: - Ready to use protocol with break-apart wells for ease of use - Standardisation and High Reproducibility - Lot to Lot Consistency -... |
| Detektionsbereich        | 0 - 1000 ng/ml* (May vary. Please confirm)                                                                                                                                                                                                                    |
| Sensitivitaet            | 30 ng/ml* (May vary. Please confirm)                                                                                                                                                                                                                          |
| Proben                   | Serum, Plasma, Cell Culture Supernatant and other biological samples. Validated for citrated/EDTA Plasma only.                                                                                                                                                |
| Anwendungsbeschreibung   | Colorimetric, Absorbance measured at 450nm                                                                                                                                                                                                                    |