

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### Recombinant Human ESAM (C-Fc) EKL-PHH0571

|                          |                                                                                                                                                                                                                                                                                    |
|--------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Artikelname              | Recombinant Human ESAM (C-Fc)                                                                                                                                                                                                                                                      |
| Artikelnummer            | EKL-PHH0571                                                                                                                                                                                                                                                                        |
| Hersteller Artikelnummer | PHH0571                                                                                                                                                                                                                                                                            |
| Alternativnummer         | EKL-PHH0571-10UG,EKL-PHH0571-50UG,EKL-PHH0571-500UG                                                                                                                                                                                                                                |
| Hersteller               | EnkiLife                                                                                                                                                                                                                                                                           |
| Wirt                     | Human                                                                                                                                                                                                                                                                              |
| Kategorie                | Proteine/Peptide                                                                                                                                                                                                                                                                   |
| Spezies Reaktivität      | Human                                                                                                                                                                                                                                                                              |
| Produktbeschreibung      | Endothelial Cell Adhesion Molecule (ESAM) is a 55 kDa type I transmembrane glycoprotein member of the JAM family of immunoglobulin superfamily molecules. The 390 amino acid Human ESAM contains a 216 amino acid extracellular domain (ECD) with a V-type...                      |
| Molekulargewicht         | 50.8 KDa                                                                                                                                                                                                                                                                           |
| NCBI                     | <a href="#">96</a>                                                                                                                                                                                                                                                                 |
| Reinheit                 | Greater than 95% as determined by reducing SDS-PAGE                                                                                                                                                                                                                                |
| Anwendungsbeschreibung   | Always centrifuge tubes before opening.Do not mix by vortex or pipetting.It is not recommended to reconstitute to a concentration less than 100µg/ml.Dissolve the lyophilized protein in distilled water.Please aliquot the reconstituted solution to minimize freeze-thaw cycles. |