

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### Recombinant CHIKV NSP2 Protein, N-GST-His, E. coli EKL-PEV60000

|                          |                                                                                                                                                                                                                                                                                    |
|--------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Artikelname              | Recombinant CHIKV NSP2 Protein, N-GST-His, E. coli                                                                                                                                                                                                                                 |
| Artikelnummer            | EKL-PEV60000                                                                                                                                                                                                                                                                       |
| Hersteller Artikelnummer | PEV60000                                                                                                                                                                                                                                                                           |
| Alternativnummer         | EKL-PEV60000-50UG,EKL-PEV60000-100UG                                                                                                                                                                                                                                               |
| Hersteller               | EnkiLife                                                                                                                                                                                                                                                                           |
| Wirt                     | E. coli                                                                                                                                                                                                                                                                            |
| Kategorie                | Proteine/Peptide                                                                                                                                                                                                                                                                   |
| Spezies Reaktivität      | Virus                                                                                                                                                                                                                                                                              |
| Produktbeschreibung      | Chikungunya virus (CHIKV) NSP2 is a multifunctional nonstructural protein essential for viral replication. It possesses RNA helicase, RNA triphosphatase, and protease activities, enabling it to unwind viral RNA, process the nonstructural polyprotein,...                      |
| Molekulargewicht         | 71.12 kDa                                                                                                                                                                                                                                                                          |
| NCBI                     | <a href="#">6</a>                                                                                                                                                                                                                                                                  |
| Reinheit                 | >90% as determined by SDS-PAGE.                                                                                                                                                                                                                                                    |
| Anwendungsbeschreibung   | Always centrifuge tubes before opening.Do not mix by vortex or pipetting.It is not recommended to reconstitute to a concentration less than 100µg/ml.Dissolve the lyophilized protein in distilled water.Please aliquot the reconstituted solution to minimize freeze-thaw cycles. |