

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Rat EndoG protein, His and GST tag, Unconjugated GTX00057-PRO

Artikelname	Rat EndoG protein, His and GST tag, Unconjugated
Artikelnummer	GTX00057-PRO
Hersteller Artikelnummer	GTX00057-pro
Alternativnummer	GTX00057-PRO-10
Hersteller	GeneTex
Kategorie	Proteine/Peptide
Applikation	FA
Spezies Reaktivität	Rat
Konjugation	Unconjugated
NCBI	362100
UniProt	Q3V5X8
Puffer	Reconstitute with 20mM Tris and 150mM NaCl to 0.1-1.0mg/ml. Do not vortex. Lyophilized from 20mM Tris, 150mM NaCl, 1mM EDTA, 1mM DTT, 0.01% SKL, 5% Trehalose, ProClin 300.
Expression System	E. coli
Formulierung	Lyophilized powder
Sequenz	Full length protein, N-terminal His and GST-Tag, Met1~Lys294 (NP_001030110.1)

Anwendungsbeschreibung	ENDOG (Endonuclease G) is a mitochondrial enzyme that cleaves DNA at double-stranded (DG)n. (DC)n and at single-stranded (DC)n tracts. This protein also has ribonuclease (RNase) and RNase H activities. It has been reported that EndoG forms complexes with some other proteins, including HSP70 (Heat shock 70 kDa protein). Besides, mouse HSPA1A (Heat shock 70 kDa protein 1A) shares similarities with rat HSP1A1 in amino acid sequence with the identity of 98.44% . Thus, a binding ELISA assay was conducted to detect the association of recombinant rat ENDOG with recombinant mouse HSPA1A. Briefly, ENDOG were diluted serially in PBS with 0.01% BSA (pH 7.4). Duplicate samples of 100 µl were then transferred to HSP1A1-coated microtiter wells and incubated for 2h at 37C. Wells were washed with PBST and incubated for 1h with anti-ENDOG pAb, then aspirated and washed 3 times. After incubation with HRP labelled secondary antibody, wells were aspirated and washed 3 times. With the addition of substrate solution, wells were incubated 15-25 minutes at 37C. Finally, add 50 µl stop solution to the wells and read at 450nm immediately. The binding activity of ENDOG with HSP1A1 was in a dose dependent manner.
------------------------	---