

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Human RPP40 protein, His tag, Unconjugated GTX00241-PRO

Artikelname	Human RPP40 protein, His tag, Unconjugated
Artikelnummer	GTX00241-PRO
Hersteller Artikelnummer	GTX00241-pro
Alternativnummer	GTX00241-PRO-10
Hersteller	GeneTex
Kategorie	Proteine/Peptide
Applikation	FA
Spezies Reaktivität	Human
Konjugation	Unconjugated
NCBI	10799
UniProt	O75818
Puffer	Reconstitute with 20mM Tris and 150mM NaCl to 0.1-1.0mg/ml. Do not vortex. Lyophilized from 20mM Tris, 150mM NaCl, 1mM EDTA, 1mM DTT, 0.01% SKL, 5% Trehalose, ProClin 300.
Expression System	E. coli
Formulierung	Lyophilized powder
Sequenz	N-terminal His-Tag, Cys49~Asn258 (NP_001273061.1)

Anwendungsbeschreibung

Ribonuclease P (RNASEP) is a type of ribonuclease which cleaves RNA. RNase P is unique from other RNases in that it is a ribozyme - a ribonucleic acid that acts as a catalyst in the same way that a protein based enzyme would. Its function is to cleave off an extra, or precursor, sequence of RNA on tRNA molecules. Besides, Methyl CpG Binding Protein 2 (MECP2) has been identified as an interactor of RNASEP, thus a binding ELISA assay was conducted to detect the interaction of recombinant human RNASEP and recombinant human MECP2. Briefly, RNASEP were diluted serially in PBS, with 0.01% BSA (pH 7.4). Duplicate samples of 100 µl were then transferred to MECP2-coated microtiter wells and incubated for 2h at 37C. Wells were washed with PBST and incubated for 1h with anti-RNASEP pAb, then aspirated and washed 3 times. After incubation with HRP labelled secondary antibody, wells were aspirated and washed 3 times. With the addition of substrate solution, wells were incubated 15-25 minutes at 37C. Finally, add 50 µl stop solution to the wells and read at 450nm immediately. The binding activity of RNASEP and MECP2 was in a dose dependent manner.