

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

SuperKine(TM) Trypsin-EDTA Solution, 0.25% (With Phenol Red) ASC-BMU109-EN

Artikelname	SuperKine(TM) Trypsin-EDTA Solution, 0.25% (With Phenol Red)
Artikelnummer	ASC-BMU109-EN
Hersteller Artikelnummer	BMU109-EN
Alternativnummer	ASC-BMU109-EN-100
Hersteller	Abbkine Scientific
Kategorie	Sonstiges
Applikation	Cell Culture
Produktbeschreibung	SuperKine(TM) Trypsin-EDTA Solution, 0.25% (With Phenol Red) contains 0.25% trypsin and 0.06% EDTA, and the pH value is 7.2-7.4. After filtration and sterilization, the digestive solution can be directly used for the digestion of cultured cells or so...
Tag	Trypsin-EDTA Solution
Formulierung	Liquid solution
Proben	Cells and Tissues
Anwendungsbeschreibung	In the process of using SuperKine(TM) Trypsin-EDTA Solution, 0.25% (With Phenol Red), special attention should be paid to avoid bacterial contamination of the digest. The digestion time of trypsin cell digestive solution should not be too long, otherwise the growth of cells will be poor after planking.