

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Human Vasodilator-stimulated phosphoprotein (VASP) ELISA Kit ASC-KTE60071

Artikelname	Human Vasodilator-stimulated phosphoprotein (VASP) ELISA Kit
Artikelnummer	ASC-KTE60071
Hersteller Artikelnummer	KTE60071
Alternativnummer	ASC-KTE60071-48, ASC-KTE60071-96
Hersteller	Abbkine Scientific
Kategorie	Kits/Assays
Applikation	ELISA
Spezies Reaktivität	Human
Produktbeschreibung	This Human Vasodilator-stimulated phosphoprotein (VASP) ELISA Kit employs a two-site sandwich ELISA to quantitate VASP....
Detektionsbereich	Please inquire
Sensitivitaet	Please inquire
Tag	VASP
NCBI	7408
UniProt	P50552
Proben	Cell culture supernatantsSerumPlasmaOther biological fluids

Anwendungsbeschreibung	This Human Vasodilator-stimulated phosphoprotein (VASP) ELISA Kit employs a two-site sandwich ELISA to quantitate VASP in samples. An antibody specific for VASP has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any VASP present is bound by the immobilized antibody. After removing any unbound substances, a biotin-conjugated antibody specific for VASP is added to the wells. After washing, Streptavidin conjugated Horseradish Peroxidase (HRP) is added to the wells. Following a wash to remove any unbound avidin-enzyme reagent, a substrate solution is added to the wells and color develops in proportion to the amount of VASP bound in the initial step. The color development is stopped and the intensity of the color is measured.
------------------------	---