

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### ELISA Flex: Bovine IFN-gamma (ALP) MAB-3119-1A-6

|                          |                                                                                                                                                                                                                                                               |
|--------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Artikelname              | ELISA Flex: Bovine IFN-gamma (ALP)                                                                                                                                                                                                                            |
| Artikelnummer            | MAB-3119-1A-6                                                                                                                                                                                                                                                 |
| Hersteller Artikelnummer | 3119-1A-6                                                                                                                                                                                                                                                     |
| Alternativnummer         | MAB-3119-1A-6                                                                                                                                                                                                                                                 |
| Hersteller               | MABTECH                                                                                                                                                                                                                                                       |
| Kategorie                | Kits/Assays                                                                                                                                                                                                                                                   |
| Applikation              | ELISA                                                                                                                                                                                                                                                         |
| Spezies Reaktivität      | Bovine, Sheep                                                                                                                                                                                                                                                 |
| Konjugation              | ALP                                                                                                                                                                                                                                                           |
| Produktbeschreibung      | ELISA Flex kit for quantification of IFN-gamma in solution, The bovine IFN-gamma ELISA Flex kit is ideal for users who want the flexibility of setting up their own ELISA assay. This product enables quantification of bovine IFN-gamma in solution, e.g.... |
| Detektionsbereich        | 4-400 pg/ml                                                                                                                                                                                                                                                   |
| NCBI                     | <a href="https://pubmed.ncbi.nlm.nih.gov/29525280">281237</a>                                                                                                                                                                                                 |
| Pubmed                   | <a href="https://pubmed.ncbi.nlm.nih.gov/29525280">29525280</a>                                                                                                                                                                                               |
| Target-Kategorie         | IFN-gamma                                                                                                                                                                                                                                                     |

|                        |                                                                                                                                                                                                                                                                                                                                           |
|------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Anwendungsbeschreibung | Complementary Products: <a href="https://www.mabtech.com/products/elisa-substrate-pnpp-alp-3652-p10-0">https://www.mabtech.com/products/elisa-substrate-pnpp-alp-3652-p10-0</a> . Contents: MT17.1,MT307, biotin,Streptavidin-ALP enzyme conjugate for ELISpot,Recombinant bovine IFN-gamma ELISA standard,Standard reconstitution buffer |
|------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|