

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

ELISpot Flex: Human Perforin (ALP) MAB-3465-2A

Artikelname	ELISpot Flex: Human Perforin (ALP)
Artikelnummer	MAB-3465-2A
Hersteller Artikelnummer	3465-2A
Alternativnummer	MAB-3465-2A
Hersteller	MABTECH
Kategorie	Kits/Assays
Applikation	ELISpot
Spezies Reaktivität	Human, Primate
Konjugation	ALP
Produktbeschreibung	Basic ELISpot kit for enumeration of cells secreting Perforin, The human Perforin ELISpot kit is ideal for users who want the flexibility of setting up their own assay. The kit enables enumeration of cells secreting human Perforin, and contains a mat...
Detektionsbereich	Streptavidin-ALP
NCBI	5551
Pubmed	25998949
Target-Kategorie	Perforin

Anwendungsbeschreibung	Complementary Products: https://www.mabtech.com/products/elispot-substrate-bcipnbt-plus-alp-3650-10-0, https://www.mabtech.com/products/anti-human-cd3-mab-cd3-2-azide-free-3605-1-0, https://www.mabtech.com/products/elispot-plates-white-3654-wp, https://www.mabtech.com/products/peppool-cef-cd8-human-3616-1-0, https://www.mabtech.com/products/peppool-cefras-global-cd8-human-3637-1-0, https://www.mabtech.com/products/peppool-cefras-global-cd4-and-cd8-human-3638-1-0. Contents: Pf-80/164, Pf-344, biotin, Streptavidin-ALP enzyme conjugate for ELISpot
------------------------	--