

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### ELISpot Plus: Human Granzyme B (ALP) MAB-3486-4APW-2

|                          |                                                                                                                                                                                                                                                               |
|--------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Artikelname              | ELISpot Plus: Human Granzyme B (ALP)                                                                                                                                                                                                                          |
| Artikelnummer            | MAB-3486-4APW-2                                                                                                                                                                                                                                               |
| Hersteller Artikelnummer | 3486-4APW-2                                                                                                                                                                                                                                                   |
| Alternativnummer         | MAB-3486-4APW-2                                                                                                                                                                                                                                               |
| Hersteller               | MABTECH                                                                                                                                                                                                                                                       |
| Kategorie                | Kits/Assays                                                                                                                                                                                                                                                   |
| Applikation              | ELISpot                                                                                                                                                                                                                                                       |
| Spezies Reaktivität      | Human                                                                                                                                                                                                                                                         |
| Konjugation              | ALP                                                                                                                                                                                                                                                           |
| Produktbeschreibung      | ELISpot Plus for enumeration of cells secreting Granzyme B, This kit is ideal for users who want a convenient and sensitive assay. The assay is designed for the enumeration of cells secreting human Granzyme B. The kit includes ELISpot plates pre-coat... |
| NCBI                     | <a href="#">3002</a>                                                                                                                                                                                                                                          |
| Target-Kategorie         | Granzyme B                                                                                                                                                                                                                                                    |

|                        |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
|------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Anwendungsbeschreibung | <p>Complementary Products:<br/> <a href="https://www.mabtech.com/products/peppool-cef-cd8-human-3616-1-0">https://www.mabtech.com/products/peppool-cef-cd8-human-3616-1-0</a>, <a href="https://www.mabtech.com/products/peppool-cefras-global-cd8-human-3637-1-0">https://www.mabtech.com/products/peppool-cefras-global-cd8-human-3637-1-0</a>, <a href="https://www.mabtech.com/products/peppool-cefras-global-cd4-and-cd8-human-3638-1-0">https://www.mabtech.com/products/peppool-cefras-global-cd4-and-cd8-human-3638-1-0</a>. Contents: MT8610, biotin,Streptavidin-ALP enzyme conjugate for ELISpot,BCIP/NBT-plus substrate for ELISpot,Pre-coated ELISpot plate, white MSIP (mAb MT28)</p> |
|------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|