

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Goat IgG anti-Rat IgG+IgM+IgA (H+L)-Biotin, MinX none DNA-SEC-183450

Artikelname	Goat IgG anti-Rat IgG+IgM+IgA (H+L)-Biotin, MinX none
Artikelnummer	DNA-SEC-183450
Hersteller Artikelnummer	SEC-183450
Alternativnummer	DNA-SEC-183450
Hersteller	dianova
Wirt	Goat
Kategorie	Antikörper
Applikation	ELISA,IHC,WB
Spezies Reaktivität	Rat
Immunogen	Anti-Rat IgG IgA IgM (H&L) was produced by repeated immunization with Rat IgG, IgA and IgM whole molecules in goat.
Konjugation	Biotin
Format	IgG
Spezifität	IgG+IgM+IgA (H+L)
Minimale Kreuzreaktivität (MinX)	no cross-adsorbtion
Produktbeschreibung	Anti-Rat IgG IgA IgM (H&L) biotin conjugated antibody generated in goat detects specifically rat IgG IgA IgM (H&L). This secondary biotin conjugated antibody anti-Rat is ideal for investigators who routinely perform titration assays, western-blot, im...

Klonalität	Polyclonal
Konzentration	1.0 mg/mL
Isotyp	Ig
Puffer	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2
Reinheit	This product was prepared from polyspecific antiserum by immunoaffinity chromatography using antigens coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities. Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-biotin, anti-Goat Serum, Rat IgA, Rat IgG and Rat IgM. This reagent is suitable for the detection of all rat immunoglobulin subclasses, isotypes and chain combinations.
Formulierung	Lyophilized
Formel	20 mM K3PO4,150 mM NaCl,pH 7,2,lyophilisate,0,01% NaN3
Target-Kategorie	Rat
Antibody Type	Secondary Antibody
Application Verdünnung	ELISA Dilution: 1:150,000, Immunohistochemistry Dilution: 1:1,000 - 1:10,000, Western Blot Dilution: 1:2,000 - 1:20,000
Anwendungsbeschreibung	Anti-Rat IgG IgA IgM (H&L) Biotin conjugated antibodies are suitable for immunoblotting (western or dot blot), ELISA, immunoperoxidase electron microscopy and immunohistochemistry as well as other peroxidase-antibody based enzymatic assays requiring lot-to-lot consistency.