

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Chicken anti Botulinum Neurotoxin Type B, Gallus NMB-AHCB

Article Name	Chicken anti Botulinum Neurotoxin Type B, Gallus
Biozol Catalog Number	NMB-AHCB
Supplier Catalog Number	AHcB
Alternative Catalog Number	NMB-AHCB
Manufacturer	NordicMubio
Host	Gallus
Category	Antikörper
Application	ELISA, WB
Species Reactivity	Bacteria
Product Description	Chicken anti Botulinum Neurotoxin Type B, Gallus 0,1 mg...
Concentration	See vial for concentration.
UniProt	P10844
Buffer	Phosphate buffered saline, pH 7.2 with 0.075% sodium azide
Source	Chickens were immunized with highly purified Botulinum Neurotoxin Type B (Clostridium botulinum). After multiple immunizations in Freund's adjuvant, eggs were collected and the IgY purified. Anti-Botulinum Neurotoxin Type B was immunoaffinity purified using
Purity	Egg-yolk derived IgY Affinity Purified
Form	Affinity purified

Formula	Phosphate buffered saline, pH 7.2 with 0.075% sodium azide
Application Notes	This antibody is suitable for ELISA and Western Blotting procedures. For direct ELISA, a 1:66,000 dilution detects Botulinum HccB and a 1:42,000 dilutions detects Clostridium botulinum Type B (BTB). For Western blot, the recommended dilution is between 1:1000 and 1:2000. Optimal working dilutions should be determined for your particular assay conditions.