

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Horse Serum Affinity Resin (For the removal of cross-reactivity to horse immunoglobulins), Unconjugated IMU-GE-096

Article Name	Horse Serum Affinity Resin (For the removal of cross-reactivity to horse immunoglobulins), Unconjugated
Biozol Catalog Number	IMU-GE-096
Supplier Catalog Number	GE-096
Alternative Catalog Number	IMU-GE-096
Manufacturer	ImmunoReagents
Category	Zellen/Zellkultur
Conjugation	Unconjugated
Concentration	>10 mg protein/ml resin
Buffer	10 mM Sodium Phosphate, 0.15 M Sodium Chloride, pH 7.2
Form	Liquid solution containing 4% agarose beads