

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

NFkB Oligonucleotide ROC-K-025

Article Name	NFkB Oligonucleotide
Biozol Catalog Number	ROC-K-025
Supplier Catalog Number	K-025
Alternative Catalog Number	ROC-K-025
Manufacturer	Rockland Immunochemicals
Category	Sonstiges
Conjugation	Unconjugated
Concentration	20 ng/μL by UV absorbance at 260 nm
Buffer	0.01 M Tris Cl, 0.15 M Sodium Chloride, 0.001 M EDTA, pH 7.4
Form	Liquid (sterile filtered)
Application Dilute	ELISA: 1:5,000 - 1:25,000, Gel Shift: See protocol, WB: 1:500 - 1:3,000
Application Notes	Gel Shift/Gel Super Shift Assays: the electrophoretic gel shift assay is used to detect sequence specific DNA-binding proteins present in nuclear extracts. For NFkB generally a HeLa nuclear extract is used. In the assay, a consensus oligonucleotide is en