

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### Colon Membrane Lysate PRS-XBL-10505

|                            |                                                                                                                    |
|----------------------------|--------------------------------------------------------------------------------------------------------------------|
| Article Name               | Colon Membrane Lysate                                                                                              |
| Biozol Catalog Number      | PRS-XBL-10505                                                                                                      |
| Supplier Catalog Number    | XBL-10505                                                                                                          |
| Alternative Catalog Number | PRS-XBL-10505-0.1                                                                                                  |
| Manufacturer               | ProSci                                                                                                             |
| Category                   | Zellen/Zellkultur                                                                                                  |
| Application                | WB                                                                                                                 |
| Concentration              | 2 mg/mL or better                                                                                                  |
| Buffer                     | HEPES pH 7.9, MgCl <sub>2</sub> , KCl, EDTA, Sucrose, Glycerol, Sodium deoxycholate, and NP-40.                    |
| Form                       | Liquid                                                                                                             |
| Application Notes          | Colon membrane protein lysate is for use in Western blotting, 10 µg to 20 µg per lane is recommended for mini gel. |