

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### GENLISA Human G Protein-Activated Inward Rectifier Potassium Channel 1 (KCNJ3 / GIRK1) ELISA KBS-KBH1846

|                            |                                                                                                                                                                                                                                                               |
|----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Article Name               | GENLISA Human G Protein-Activated Inward Rectifier Potassium Channel 1 (KCNJ3 / GIRK1) ELISA                                                                                                                                                                  |
| Biozol Catalog Number      | KBS-KBH1846                                                                                                                                                                                                                                                   |
| Supplier Catalog Number    | KBH1846                                                                                                                                                                                                                                                       |
| Alternative Catalog Number | KBS-KBH1846-1                                                                                                                                                                                                                                                 |
| Manufacturer               | Krishgen Biosystems                                                                                                                                                                                                                                           |
| Category                   | Kits/Assays                                                                                                                                                                                                                                                   |
| Species Reactivity         | Human                                                                                                                                                                                                                                                         |
| Product Description        | The GENLISA Human G Protein-Activated Inward Rectifier Potassium Channel 1 (KCNJ3 / GIRK1) ELISA includes features like: - Ready to use protocol with break-apart wells for ease of use - Standardisation and High Reproducibility - Lot to Lot Consistenc... |
| Range                      | 0 - 1000 ng/ml* (May vary. Please confirm)                                                                                                                                                                                                                    |
| Sensitivity                | 30 ng/ml* (May vary. Please confirm)                                                                                                                                                                                                                          |
| Samples                    | Serum, Plasma, Cell Culture Supernatant and other biological samples. Validated for citrated/EDTA Plasma only.                                                                                                                                                |
| Application Notes          | Colorimetric, Absorbance measured at 450nm                                                                                                                                                                                                                    |