

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### Mouse TIMP1 (Tissue Inhibitors Of Metalloproteinase 1) ELISA Kit EKL-EM22719

|                            |                                                                                                                                                                                                                                                               |
|----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Article Name               | Mouse TIMP1 (Tissue Inhibitors Of Metalloproteinase 1) ELISA Kit                                                                                                                                                                                              |
| Biozol Catalog Number      | EKL-EM22719                                                                                                                                                                                                                                                   |
| Supplier Catalog Number    | EM22719                                                                                                                                                                                                                                                       |
| Alternative Catalog Number | EKL-EM22719-96T                                                                                                                                                                                                                                               |
| Manufacturer               | EnkiLife                                                                                                                                                                                                                                                      |
| Category                   | Kits/Assays                                                                                                                                                                                                                                                   |
| Species Reactivity         | Mouse                                                                                                                                                                                                                                                         |
| Product Description        | This ELISA kit uses the Sandwich-ELISA principle. The micro ELISA plate provided in this kit has been pre-coated with an antibody specific to Mouse TIMP1. Standards or samples are added to the micro ELISA plate wells and combined with the specific an... |
| Range                      | 78-5000pg/mL                                                                                                                                                                                                                                                  |
| Sensitivity                | 28pg/mL                                                                                                                                                                                                                                                       |
| Samples                    | Serum,plasma and other biological fluids                                                                                                                                                                                                                      |
| Target                     | TIMP1                                                                                                                                                                                                                                                         |