

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### Porcine IP-10/CXCL10 (Interferon Gamma Induced Protein 10kDa) ELISA Kit EKL-EP27621

|                            |                                                                                                                                                                                                                                                               |
|----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Article Name               | Porcine IP-10/CXCL10 (Interferon Gamma Induced Protein 10kDa) ELISA Kit                                                                                                                                                                                       |
| Biozol Catalog Number      | EKL-EP27621                                                                                                                                                                                                                                                   |
| Supplier Catalog Number    | EP27621                                                                                                                                                                                                                                                       |
| Alternative Catalog Number | EKL-EP27621-96T                                                                                                                                                                                                                                               |
| Manufacturer               | EnkiLife                                                                                                                                                                                                                                                      |
| Category                   | Kits/Assays                                                                                                                                                                                                                                                   |
| Species Reactivity         | Porcine                                                                                                                                                                                                                                                       |
| Product Description        | This ELISA kit uses the Sandwich-ELISA principle. The micro ELISA plate provided in this kit has been pre-coated with an antibody specific to Porcine IP-10/CXCL10. Standards or samples are added to the micro ELISA plate wells and combined with the sp... |
| Range                      | 31.25-2000pg/mL                                                                                                                                                                                                                                               |
| Sensitivity                | 18.75pg/mL                                                                                                                                                                                                                                                    |
| Samples                    | Serum,plasma and other biological fluids                                                                                                                                                                                                                      |
| Target                     | IP-10/CXCL10                                                                                                                                                                                                                                                  |