

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### **GENLISA Human Phospho-C-Jun N-Terminal Kinases / stress-Activated Protein Kinase (P-JNK / SAPK) ELISA KBS-KBH11252**

Article Name	GENLISA Human Phospho-C-Jun N-Terminal Kinases / stress-Activated Protein Kinase (P-JNK / SAPK) ELISA
Biozol Catalog Number	KBS-KBH11252
Supplier Catalog Number	KBH11252
Alternative Catalog Number	KBS-KBH11252-1
Manufacturer	Krishgen Biosystems
Category	Kits/Assays
Species Reactivity	Human
Product Description	The GENLISA Human Phospho-C-Jun N-Terminal Kinases / stress-Activated Protein Kinase (P-JNK / SAPK) ELISA includes features like: Ready to use protocol with break-apart wells for ease of use Standardisation and High Reproducibility Lot to Lot Consist...
Range	0 - 1000 ng/ml* (May vary. Please confirm)
Sensitivity	30 ng/ml* (May vary. Please confirm)
Samples	Serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids. For tear samples, please connect with us to discuss further.
Application Notes	Colorimetric, Absorbance measured at 450nm