

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### Human FBLN5 (Fibulin 5) ELISA Kit EKL-EH24962

|                            |                                                                                                                                                                                                                                                               |
|----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Article Name               | Human FBLN5 (Fibulin 5) ELISA Kit                                                                                                                                                                                                                             |
| Biozol Catalog Number      | EKL-EH24962                                                                                                                                                                                                                                                   |
| Supplier Catalog Number    | EH24962                                                                                                                                                                                                                                                       |
| Alternative Catalog Number | EKL-EH24962-96T                                                                                                                                                                                                                                               |
| Manufacturer               | EnkiLife                                                                                                                                                                                                                                                      |
| Category                   | Kits/Assays                                                                                                                                                                                                                                                   |
| Species Reactivity         | Human                                                                                                                                                                                                                                                         |
| Product Description        | This ELISA kit uses the Sandwich-ELISA principle. The micro ELISA plate provided in this kit has been pre-coated with an antibody specific to Human FBLN5. Standards or samples are added to the micro ELISA plate wells and combined with the specific an... |
| Range                      | 15.6-1000pg/mL                                                                                                                                                                                                                                                |
| Sensitivity                | 5.6pg/mL                                                                                                                                                                                                                                                      |
| Samples                    | Serum,plasma and other biological fluids                                                                                                                                                                                                                      |
| Target                     | FBLN5                                                                                                                                                                                                                                                         |