

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### Rat TAP (Trypsinogen Activation Peptide) ELISA Kit EKL-ER21443

|                            |                                                                                                                                                                                                                                                               |
|----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Article Name               | Rat TAP (Trypsinogen Activation Peptide) ELISA Kit                                                                                                                                                                                                            |
| Biozol Catalog Number      | EKL-ER21443                                                                                                                                                                                                                                                   |
| Supplier Catalog Number    | ER21443                                                                                                                                                                                                                                                       |
| Alternative Catalog Number | EKL-ER21443-96T                                                                                                                                                                                                                                               |
| Manufacturer               | EnkiLife                                                                                                                                                                                                                                                      |
| Category                   | Kits/Assays                                                                                                                                                                                                                                                   |
| Species Reactivity         | Rat                                                                                                                                                                                                                                                           |
| Product Description        | This ELISA kit uses the Competitive-ELISA principle. The micro ELISA plate provided in this kit has been pre-coated with Rat TAP. During the reaction, Rat TAP in the sample or standard competes with a fixed amount of Rat TAP on the solid phase suppor... |
| Range                      | 123.5-10000pg/mL                                                                                                                                                                                                                                              |
| Sensitivity                | 50.3pg/mL                                                                                                                                                                                                                                                     |
| Samples                    | Serum,plasma and other biological fluids                                                                                                                                                                                                                      |
| Target                     | TAP                                                                                                                                                                                                                                                           |