

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Human SEMA3F protein, His tag (active), Unconjugated GTX00267-PRO

Article Name	Human SEMA3F protein, His tag (active), Unconjugated
Biozol Catalog Number	GTX00267-PRO
Supplier Catalog Number	GTX00267-pro
Alternative Catalog Number	GTX00267-PRO-10
Manufacturer	GeneTex
Category	Proteine/Peptide
Application	FA
Species Reactivity	Human
Conjugation	Unconjugated
NCBI	6405
UniProt	Q13275
Buffer	Reconstitute with 20mM Tris and 150mM NaCl to 0.1-1.0mg/ml. Do not vortex. Lyophilized from 20mM Tris, 150mM NaCl, 1mM EDTA, 1mM DTT, 0.01% SKL, 5% Trehalose, ProClin 300.
Expression System	E. coli
Form	Lyophilized powder
Sequence	N-terminal His-Tag, Val350~Leu483 (NP_001305729.1)

Application Notes	Semaphorin-3F (SEMA3F) is a member of semaphorins family. All the family members have a secretion signal, a 500-amino acid sema domain, and 16 conserved cysteine residues. Semaphorins not only guide axons in development, but also have major roles in immune function (classes 4, 6, and 7) and the development of bones. Class 3 semaphorins are one of the most versatile semaphorin classes, in which Sema3a is the most studied. Besides, SEMA3F is deleted or downregulated in many metastatic tumors. To test the effect of SEMA3F on inhibit cell proliferation, human lung carcinoma cells A549 were seeded into triplicate wells of 96-well plates at a density of 5000 cells/well and allowed to attach, replaced with serum-free overnight, then the medium was replaced with 2% serum standard DMEM prior to the addition of various concentrations of recombinant human SEMA3F. After incubated for 48h, cells were observed by inverted microscope and cell proliferation was measured by Cell Counting Kit-8 (CCK-8). Briefly, 10 μl of CCK-8 solution was added to each well of the plate, then the absorbance at 450nm was measured using a microplate reader after incubating the plate for 1-4 hours at 37C. Proliferation of A549 cells after incubation with SEMA3F for 48h observed by inverted microscope. Cell viability was assessed by CCK-8 assay after incubation with recombinant SEMA3F for 48h. And SEMA3F significantly increased cell viability of A549 cells.
-------------------	---