

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### High Sensitive ELISA Kit for Immunoglobulin G (IgG) USC-HEA544MU

|                            |                                                                                                      |
|----------------------------|------------------------------------------------------------------------------------------------------|
| Article Name               | High Sensitive ELISA Kit for Immunoglobulin G (IgG)                                                  |
| Biozol Catalog Number      | USC-HEA544MU                                                                                         |
| Supplier Catalog Number    | HEA544Mu                                                                                             |
| Alternative Catalog Number | USC-HEA544MU-48, USC-HEA544MU-96                                                                     |
| Manufacturer               | Cloud-Clone                                                                                          |
| Category                   | Kits/Assays                                                                                          |
| Species Reactivity         | Mouse                                                                                                |
| Product Description        | 48T, 96T, 96T*5, 96T*10, 96T*100...                                                                  |
| Range                      | 0.78-50ng/mL                                                                                         |
| Sensitivity                | The minimum detectable dose of this kit is typically less than 0.26ng/mL                             |
| Samples                    | serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids |
| Application Notes          | Enzyme-linked immunosorbent assay for Antigen Detection.                                             |