

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### **Rabbit Antiserum anti-Dog IgG (F(ab)2)-unconj., MinX none DNA-SEC-182420**

|                            |                                                           |
|----------------------------|-----------------------------------------------------------|
| Article Name               | Rabbit Antiserum anti-Dog IgG (F(ab)2)-unconj., MinX none |
| Biozol Catalog Number      | DNA-SEC-182420                                            |
| Supplier Catalog Number    | SEC-182420                                                |
| Alternative Catalog Number | DNA-SEC-182420                                            |
| Manufacturer               | dianova                                                   |
| Host                       | Rabbit                                                    |
| Category                   | Antikörper                                                |
| Application                | ELISA,IHC,WB                                              |
| Species Reactivity         | Canine                                                    |
| Immunogen                  | Dog IgG F(ab)2 fragment, rabbit anti Dog Fab2             |
| Conjugation                | Unconjugated                                              |
| Format                     | Antiserum                                                 |
| Target Specificity         | IgG (F(ab')2)                                             |
| Cross-Adsorption (MinX)    | no cross-adsorbtion                                       |
| Clonality                  | Polyclonal                                                |
| Concentration              | 32.0 mg/mL                                                |
| Buffer                     | 0.01 M Sodium Phosphate, 0.15 M Sodium Chloride, pH 7.2   |

|                    |                                                                                                                                                                                                                                                                                             |
|--------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Purity             | This product was prepared from monospecific antiserum by a delipidation and defibrination. Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-rabbit serum, Dog IgG, Dog IgG F(ab') <sub>2</sub> and Dog Serum. No reaction was observed against Dog IgG F(c). |
| Form               | Lyophilized                                                                                                                                                                                                                                                                                 |
| Antibody Type      | Secondary Antibody                                                                                                                                                                                                                                                                          |
| Application Dilute | ELISA Dilution: 1:20,000 - 1:100,000, Immunohistochemistry Dilution: 1:1,000 - 1:5,000, Western Blot Dilution: 1:2,000 - 1:10,000                                                                                                                                                           |
| Application Notes  | Secondary antibody reagents are ideal for ELISA, western blotting, Immunohistochemistry, Fluorescence Microscopy, Flow Cytometry as well as other antibody detection methods.                                                                                                               |