

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### **Ab-10 Rapid HiLyte Fluor 647 Labeling Kit DOJ-LK36-10**

|                            |                                                                                                                                                                                       |
|----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Article Name               | Ab-10 Rapid HiLyte Fluor 647 Labeling Kit                                                                                                                                             |
| Biozol Catalog Number      | DOJ-LK36-10                                                                                                                                                                           |
| Supplier Catalog Number    | LK36-10                                                                                                                                                                               |
| Alternative Catalog Number | DOJ-LK36-10                                                                                                                                                                           |
| Manufacturer               | Dojindo                                                                                                                                                                               |
| Category                   | Kits/Assays                                                                                                                                                                           |
| Product Description        | - Small Quantity of Antibody Requirement : 10 ug Antibody Sample Labeling - Easy Labeling Procedure : Just Mix Antibody and Labeling Agent - Rapid Labeling : Less than 30 minutes... |