

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Tissue, Section, Human Disease, Arteriosclerosis, Heart, Aorta and Matched Non-arteriosclerosis Aorta (Frozen) HisTek(TM) USB-T5595-5112G2

| | |
|----------------------------|---|
| Article Name | Tissue, Section, Human Disease, Arteriosclerosis, Heart, Aorta and Matched Non-arteriosclerosis Aorta (Frozen) HisTek(TM) |
| Biozol Catalog Number | USB-T5595-5112G2 |
| Supplier Catalog Number | T5595-5112G2 |
| Alternative Catalog Number | USB-T5595-5112G2-2X5 |
| Manufacturer | US Biological |
| Category | Sonstiges |
| Product Description | Tissues were snap frozen in liquid nitrogen immediately after excision, and embedded in OCT. All frozen tissue products, including frozen tissue arrays, frozen tissue sections, and frozen tissue panels, can be used for both immunohistochemistry and i... |