

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

mIHC Complete TSA Multiplex Kit - 10 Color Panel ABC-A334673-100

Article Name	mIHC Complete TSA Multiplex Kit - 10 Color Panel
Biozol Catalog Number	ABC-A334673-100
Supplier Catalog Number	A334673-100
Alternative Catalog Number	ABC-A334673-100
Manufacturer	Antibodies.com
Category	Kits/Assays
Application	IHC
Product Description	Multiplex Fluorescence Immunohistochemistry (mIHC) Kit with TSA Amplification with 10 colors, full-spectrum panel [440, 470, 520, 540, 570, 620, 650, 660, 670, 700]....
Form	Liquid. The Repair Quick Removal 2-in-1 Buffer in this kit is powder.
Samples	This kit is suitable for paraffin-embedded tissues, frozen tissues, tissue microarrays (TMA), and cell climbing slides.
Application Dilute	The fluorescent dye is provided as a 100* stock solution. The working solution is prepared by diluting the stock 1:100 in the freshly prepared signal amplification reaction solution. TSA470SN, a nuclear staining dye, should be diluted 1:100 in sterile wat