

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### **2 PRO Screen, 2 Pressure sealing plate, 4 wash/eluate collection plates (2 ml), 4 eluate collection plates (0.35 ml), 4 blotting membranes, 2 sealing mats** **PAR-PAL-HT-Q-2**

|                            |                                                                                                                                                               |
|----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Article Name               | 2 PRO Screen, 2 Pressure sealing plate, 4 wash/eluate collection plates (2 ml), 4 eluate collection plates (0.35 ml), 4 blotting membranes, 2 sealing mats    |
| Biozol Catalog Number      | PAR-PAL-HT-Q-2                                                                                                                                                |
| Supplier Catalog Number    | PAL-HT-Q-2                                                                                                                                                    |
| Alternative Catalog Number | PAR-PAL-HT-Q-2                                                                                                                                                |
| Manufacturer               | ProteinArk                                                                                                                                                    |
| Category                   | Kits/Assays                                                                                                                                                   |
| Product Description        | 2 PRO Screen, 2 Pressure sealing plate, 4 wash/eluate collection plates (2 ml), 4 eluate collection plates (0.35 ml), 4 blotting membranes, 2 sealing mats... |