

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### ELISA Flex: Porcine IFN-gamma (ALP) MAB-3130-1A-6

|                            |                                                                                                                                                                                                                                                               |
|----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Article Name               | ELISA Flex: Porcine IFN-gamma (ALP)                                                                                                                                                                                                                           |
| Biozol Catalog Number      | MAB-3130-1A-6                                                                                                                                                                                                                                                 |
| Supplier Catalog Number    | 3130-1A-6                                                                                                                                                                                                                                                     |
| Alternative Catalog Number | MAB-3130-1A-6                                                                                                                                                                                                                                                 |
| Manufacturer               | MABTECH                                                                                                                                                                                                                                                       |
| Category                   | Kits/Assays                                                                                                                                                                                                                                                   |
| Application                | ELISA                                                                                                                                                                                                                                                         |
| Species Reactivity         | Porcine                                                                                                                                                                                                                                                       |
| Conjugation                | ALP                                                                                                                                                                                                                                                           |
| Product Description        | ELISA Flex kit for quantification of IFN-gamma in solution, The porcine IFN-gamma ELISA Flex kit is ideal for users who want the flexibility of setting up their own ELISA assay. This product enables quantification of porcine IFN-gamma in solution, e.... |
| Range                      | 16-1600 pg/ml                                                                                                                                                                                                                                                 |
| NCBI                       | <a href="#">396991</a>                                                                                                                                                                                                                                        |
| Pubmed                     | <a href="https://pubmed.ncbi.nlm.nih.gov/25036760">https://pubmed.ncbi.nlm.nih.gov/25036760</a> target=_blank rel=nofollow>25036760                                                                                                                           |
| Target                     | IFN-gamma                                                                                                                                                                                                                                                     |

|                   |                                                                                                                                                                                                                                                                                                                                                  |
|-------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Application Notes | Complementary Products: <a href="https://www.mabtech.com/products/elisa-substrate-pnpp-alp-3652-p10-0">https://www.mabtech.com/products/elisa-substrate-pnpp-alp-3652-p10-0</a> . Contents: pIFNgamma-I,P2C11, biotin ,Streptavidin-ALP enzyme conjugate for ELISpot,Recombinant porcine IFN-gamma ELISA standard,Standard reconstitution buffer |
|-------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|